Test na kotinin, nikotin (ze slin)

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REF DCT-803	Čeština
-------------	---------

[SHRnutí]

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Koncentrace (Cut-off Range)	n	COT	
		-	+
0% Cut-off	30	30	0
-50% Cut-off	30	30	0
-25% Cut-off	30	25	5
Cut-off	30	20	10
+25% Cut-off	30	7	23
+50% Cut-off	30	0	30
+300% Cut-off	30	0	30

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Amobarbital	Efavirenz (Sustiva)	l-Methamphetamine	Prednisolone
Amoxapine	EMDP	Methaqualone	Prednisone
Amoxicillin	Ephedrine	Methoxyphenamine	5β-Pregnane-3α,17α,21α-triol
l-Amphetamine	l-Ephedrine	(-) 3,4-Methylenedioxyamphetamine (MDA)	Procaine
Ampicillin	(±)-Epinephrine	(+) 3,4 Methylenedioxy-methamphetamine (MDMA)	Promazine
Apomorphine	l-Epinephrine		Promethazine
l-Ascorbic acid	Erythromycin		d,l-Propanolol
Aspartame	β-Estradiol		d-Propoxyphene
Atropine	Estrone-3-sulfate	Methylphenidate	d-Pseudoephedrine
Benzilic acid	Ethanol (Ethyl alcohol)	Methpyrion	Quinacrine
Benzoic acid	Ethyl-p-aminobenzoate	Methaqualone	Quinidine
Benzoylcegonine	Etodolac	Metoprolol	Quinine
Benztphetamine	Famprophazone	Morphine sulfate	Ranitidine
Bilirubin	Fenfluramine	Morphine-	Riboflavin
(±)-Brompheniramine	Fenoprofen	3-β-D-glucuronide	Salicylic acid
Buspirone	Fentanyl	Nalidixic acid	Secobarbital
Caffeine	Fluoxetine	Nalorphine	Serotonin
Cannabidiol	Furosemide	Naloxone	(5-hydroxytryptamine)
Cannabinol	Genitistic acid	Naltrexone	Sodium chloride
Chloral hydrate	d (+) Glucose	Methpyrion	Sulfamethazine
Chloramphenicol	Guaiacol glyceryl ether	Metoprolol	Sulindac
Chlorphenazexide	Guaiacol glyceryl ether carbamate	Nimesulide	Temazepam
Chloroquine	Hemoglobin	Norcodein	Tetracycline
Chlorothiazide		Morphine sulfate	Tetrahydrocortisone,
(+)-Chlorpheniramine	Hydralazine	α-Naphthaleneacetic acid	3-acetate
(±)-Chlorpheniramine	Hydrochlorothiazide	Norethindrone	Tetrahydrozoline
Chlorpromazine	Hydrocodone	Thiobaine	Thebaine
Chlorprothixene	Hydrocortisone	d-Norpropoxyphene	Theophylline
Cholesterol	Hydromorphone	Noscapine	Thiamine
Cimetidine	p-Hydroxyamphetamine	d,l-Octopamine	Thioridazine
Clomipramine	o-Hydroxyhippuric acid	Orphenadrine	(chlorpromazine)
Clonidine	p-Hydroxymethamphetamine	Oxalic acid	l-Thyroxine
Cocaine	p-Hydroxynorephedrine	Oxazepam	Tolbutamide
Codeine	Hydroxyzine	Oxolinic acid	cis-Tramadol
Cortisone	3-Hydroxytyramine	Oxycodone	Trazodone
Creatinine	Ibuprofen	Oxymetazoline	Triamterene
Cyclobarbitol	Imipramine	Oxymorphone	Trifluoperazine
Cyclobenzaprine	Iproniazid	Papaverine	Trimethobenzamide
Deoxycorticosterone	(-) Isoproterenol	Pemoline	Trimethoprim
(-) Deoxyephedrine	Isoxsuprine	Penicillin-G	Trimipramine
R (-) Deprenyl	Kanamycin	Pentazocine	Tryptamine
Dextromethorphan	Ketamine	Pentobarbital	d,l-Tryptophan
Diazepam	Ketoprofen	Perphenazine	Tyramine
Diclofenac	Labetalol	Phencyclidine	d,l-Tyrosine
Dicyclomine	Levorphanol	Phenelzine	Uric acid
Diflunisal	Lidocaine	Pheniramine	Verapamil
Digoxin	Lindane	Phenobarbital	Zomepirac
	(hexachlorocyclohexane)	Phenothiazine	

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